

Cephas Medical's advanced Amnicot design helps women avoid instrumental rupture of the membrane. The Amnicot is a single use fingertip with a plastic hook for performing an artificial rupture of the amniotic membrane (amniotomy)

What is it used for?

It is used for amniotomy, also known as artificial rupture of membranes (AROM) and by the lay description "breaking the water". This is the intentional rupture of the amniotic sac by an obstetrical provider.



PARAMETERS	SPECIFICATIONS
Material	Nitrile & Natural rubber
Hook	Medical grade Polycarbonate
Sterilization	EO sterile
Size	One size fits for all

Features

- One finger usage and no other instrument is inserted through the cervix
- The use of instruments can cause significant anxiety to patients
- It is possible to feel exactly where the hook is applied
- The doctor can stabilize the presenting part with the other hand
- As finger is flexible, instrument lithotomy is not necessary



Less Pain



Less Anxiety



Exact Placement
of Hooks



No Assistance
Required



No Lithotomy
Necessary



USFDA Registration
Certificate



ISO 13485 : 2016



Facility: Special Economic
Zone Chennai, India.

Manufacturer :
Cephas Medical Private Limited

📍 B13, MEPZ Special Economic Zone,
Chennai - 600 045, India.

✉ cephas@cephasmedical.net

☎ +91 98403 30964
+91 87544 73188